



Hearing aid management

To ensure you receive maximum benefit from your hearing device, we recommend following these care & maintenance and usage tips.

1. Leave your hearing aids in all day. Only take them out when:

- Going to bed
- Going for a shower or bath

2. In the morning:

- Make sure the battery is in the hearing aids
- Close the battery door
- Put hearing aids in
- **Red** goes in the **right** ear
- **Blue** goes in the **left** ear

3. At night:

- Take hearing aids out before bed
- Open the battery doors
- Place hearing aid in a case
- No need to shut the case

4. Weekly:

- Change the battery on a specified day each week
- Clean and or replace the wax busters weekly (more if it gets blocked)
- Brush over the microphone covers



Putting the batteries in the hearing aids:

1. Open the battery door
2. Take the a battery out of the packed and take the sticker off
3. Put the battery in the hearing aid with the flat side facing up

If the hearing aid is dead/does not want to work:

1. Change the battery
2. Check the wax buster
3. Brush over the microphones and where the sound comes out
4. If it still does not work – please bring the hearing aids in to us.